

## TECHNICAL FEATURES

- MODEL: FISIOFIELD MAXI.
- COMMERCIAL CLASSIFICATION: magnetotherapy device.
- TECHNICAL CLASSIFICATION: electromedical device, Class I type BF.
- MEDICAL DEVICE CLASS: IIa (Dir. 93/42/EEC, modified by the Directive 2007/47/EC).
- REFERENCE STANDARDS: EN 60601-1 (IEC 60601-1), EN 60601-1-2 (IEC 60601-1-2).
- GENERATOR SIZE: 320x430x980 mm.
- WEIGHT: 28 kg.
- POWER SUPPLY VOLTAGE: 230V single phase (115V upon request).
- LINE FREQUENCY: 50-60Hz.
- ABSORBED POWER: 400VA.
- WAVE SHAPES: square and sinusoidal, selectable.
- EMISSION FREQUENCY:
  - Square wave: 5-100Hz (E.L.F.) and 100-1000Hz (Low and Medium frequency).
  - Sinusoidal wave: 1-100 Hz.
- PEAK INTENSITY OF THE MAGNETIC FIELD:
  - $\geq 100$  Gauss each applicator.
  - $\geq 200$  Gauss each pair of contact applicators  $\varnothing$  9 cm and  $\varnothing$  13 cm.
- MEDIUM INTENSITY OF THE MAGNETIC FIELD:
  - $\geq 50$  Gauss each applicator, adjustable from 0 to 100% in increments of 10.
  - $\geq 100$  Gauss each pair of contact applicators  $\varnothing$  9 cm and  $\varnothing$  13 cm, adjustable from 0 to 100% in increments of 10.
- NUMBER OF OUTPUTS: 4 totally independent.
- DISPLAY:
  - Backlit graphic display
  - 7" Color Touch Screen Display OPTIONAL
- TIMER: ranging from 1 min up to 95 min.
- ACOUSTIC SIGNAL at end of treatment.
- PRESET THERAPEUTIC PROGRAMS: 50 for the main pathologies.
- USER'S PROGRAMS: customizable programs.
- STANDARD ACCESSORIES:
  - Therapy bed with solenoid  $\varnothing$  60 cm.
  - Test magnet.
  - Power supply cable.
- OPTIONAL ACCESSORIES:
  - Additional therapy bed with solenoid  $\varnothing$  60 cm or  $\varnothing$  80 cm.
  - Additional solenoid  $\varnothing$  60 cm or  $\varnothing$  80 cm for bed (Total body).
  - Portable solenoid for limbs  $\varnothing$  30 cm.
  - Pair of contact applicators  $\varnothing$  9 cm.
  - Pair of contact applicators  $\varnothing$  13 cm.
  - Mats for magnetotherapy of different sizes.
  - Flexible applicators.
  - Elastic straps.
- THERAPY BED dimensions: 198x65x75 cm.
- **CE** MARK : device complying with Directive 93/42/EEC, modified by the Directive 2007/47/EC.

 FISIOLINE® instruments are covered by UNIPOLSAI product liability insurance.



MADE IN ITALY

### Fisioline® srl

Borgata Molino, 29 • 12060 VERDUNO (CN) • ITALY

Tel.: +39.0172.470432-0172.470433 • Fax.: +39.0172.470891

<http://www.fisioline.com> • e-mail: [fisioline@fisioline.com](mailto:fisioline@fisioline.com)

Technical features and sizes of the devices described herein may be subject to change, based on technological evolution and constant improvements to production.



Ed. 10/20



# FISIOFIELD MAXI

MAGNETOTHERAPY

